

Biological Indicators

Steam

- Self-contained Biological Indicators
- Readout: 24Hours.
- Geobacillus Stearothermophilus (#7953), 10⁶ Spores per vial.
- Reliable routine monitoring.
- Identification cap color.



Biological Indicators

Hydrogen Peroxide

- Self-contained Biological Indicators
- Readout: 24Hours.
- Geobacillus Stearothermophilus (#7953), 10⁶ Spores per vial.
- Reliable routine monitoring.
- Identification cap color.



Chemical Indicators

Type 4 Indicator

- Steam multiparameter indicator.
- Lead free
- 160 x 15 mm strips, Bag per 250 units



Chemical Indicators

Bowie&Dick Test Sheets

- Bowie and Dick test sheet for verification of steam penetration
- Lead free.
- Final turning to black.
- Box per 50 sheets.



Chemical Indicators

Disposable Bowie&Dick Test pack

- Uniform, barrier in each pack.
- Lead free.
- Final turning to black.
- Pack of 120x120x15mm, Box per 20 units.

