



- Switching ON/OFF
- Illumination intensity
- Color temperature adjustment
- Endo-mode
- Standard illumination mode
- Sofer light mode
- Synchronic mode(option)
- Camera model(option)



Using time of battery: 2.5 hours
Battery charging time: 4 hours



Technical data

Diameter of lamphead:	500mm
Numbers of LEDs:	48LEDs
LED life span:	50000hrs
Power supply:	AC 110~220V 50/60Hz
Total power input:	65VA
Central illuminance(EC):	40kLx-160kLx
Brightness control:	0%-100%
Light field diameter (D10) :	120-280mm
Light spot distribution(D50 / D10):	0.5
Depth of illumination:	600mm
Color teperature:	3500-5000K
Color rendering index:	100>Ra>85
Height of instanllation:	250-320 cm

Battery system running automatically
in case of electricity failure and providing
continuous energy for >3 hours

CE

